



# General Chemistry Laboratory

---

## Quantitative Analysis of Cobalt(II) Ions



# Preparation

## Collect the following items

- Six 50 mL Erlenmeyer flasks + Six cork stoppers
- One 10 mL volumetric flask
- One 2 mL graduated pipet and one pipet filler
- One cuvette + Two lens tissues
- One drop pipet



## From your personal equipment

- Use one clean and dry test tube to take 8 mL standard  $\text{Co}^{2+}$  solution
- 10 mL and 50 mL graduated cylinders
- The TA will distribute a  $\text{Co}^{2+}$  solution with unknown concentration to each group



# Objective and Principles

- **Objective:** Determine trace amounts of  $\text{Co}^{2+}$  ions by quantitative analysis of  $[\text{Co}(\text{SCN})_4]^{2-}$  with a spectrophotometer





- **Lab techniques:**
  - Operate a spectrophotometer
  - Volumetric flask
  - Graduated pipet



# Color and Absorption of $[\text{Co}(\text{H}_2\text{O})_6]^{2+}$

- $\text{Co}^{2+}$  has 7 electrons in the 3d orbitals
- In aqueous environment, each  $\text{Co}^{2+}$  is coordinated by six  $\text{H}_2\text{O}$  ligands (oxygen atom donates lone pair electrons) and forms a  $[\text{Co}(\text{H}_2\text{O})_6]^{2+}$  complex ion
- The energy of 3d orbitals are splitted in the octahedral crystal field
- Low absorption probability

A compound's visible color is complementary to the color it absorbs:







# Color and Absorption of $[\text{Co}(\text{SCN})_4]^{2-}$

- Each  $\text{Co}^{2+}$  is coordinated by four  $\text{SCN}^-$  ligands (nitrogen atom donates lone pair electrons) and forms a tetrahedral  $[\text{Co}(\text{SCN})_4]^{2-}$  complex ion
- The energy of 3d orbitals are splitted differently in this tetrahedral crystal field (2 orbitals have lower energy than the other three orbitals)
- Absorbs red/orange light with higher probability (stronger absorption)

A compound's visible color is complementary to the color it absorbs:






# Beer's Law


$$A(\lambda) = \epsilon(\lambda) \times b \times c$$

A: Measured absorbance  
 $\epsilon$ : Extinction coefficient  
b: Path length  
c: Concentration of the absorbing substance



- Transmittance:  $T = P/P_0$
- Absorbance:  $A = -\log T$   
 $A = -\log(P/P_0)$

- For a dilute solution, the absorbance is proportional to the concentration
- The extinction coefficient  $\epsilon$  can be calculated from the Abs vs. concentration calibration curve
- With  $\epsilon$  known, the unspecified conc. of sample solution can be determined from its absorbance





# Step 1: Calibration of Spectrophotometer

- 1) Switch on the power and let warm up for at least 20 minutes
- 2) Check whether the cuvette holder is empty, then close the lid of sample compartment
- 3) Press the “**A/T/C**” button to set the measurement mode to A (absorbance)
- 4) Set the analytical wavelength to 620 nm
- 5) Press the “**BLANK**” button to zero the absorbance reading

**(4) Wavelength setting**

**(1) Power switch**

**(3) A/T/C Mode**

**(5) Blank**





**(2) Cuvette holder**



# Step 2: Prepare Blank Solution

- Measure 0.8 mL 6 M HCl, 2.0 mL 50% KSCN, and 4.8 mL acetone into a clean 10 mL volumetric flask
- Add DI water until the liquid level matches the inscribed mark; Install the stopper cap (hold with a finger), invert the flask several times to ensure thorough mixing
- Transfer the solution to a clean and dry 50 mL Erlenmeyer flask, and seal it with a cork stopper (acetone is a volatile solvent)



✓ White precipitates may appear depending on the sequence of mixing chemicals (why?)





# Step 3: Prepare $\text{Co}^{2+}$ Solutions

- According to the values listed in the table below, use a 2 mL graduated pipet to accurately measure 0.50 – 2.00 mL standard  $\text{Co}^{2+}$  solution (0.10 mg/mL) and transfer it to the clean 10 mL volumetric flask
- Add 0.8 mL 6 M HCl, 2.0 mL 50% KSCN, and 4.8 mL acetone into the same volumetric flask
- Dilute each solution to 10 mL with DI water, mix thoroughly, and transfer to clean and dry Erlenmeyer flasks (sealed with cork stoppers)



	Sample No.	0.10 mg/mL $\text{Co}^{2+}$ (mL)	6 M HCl	50% KSCN	Acetone
Blank solution (step 2) →	0	0	0.8 mL	2.0 mL	4.8 mL
For calibration curve {	1	0.50			
	2	1.00			
	3	1.50			
	4	2.00			

✓ Use the same volumetric flask to prepare all solutions in the sequence 0 → 4




## Step 4: Prepare Unknown $\text{Co}^{2+}$ Sample

- Use a 2 mL graduated pipet to transfer  $x$  mL of the unknown conc.  $\text{Co}^{2+}$  solution (assigned to your group) into a clean 10 mL volumetric flask
- Add 0.8 mL 6 M HCl, 2.0 mL 50% KSCN, and 4.8 mL acetone into the same volumetric flask
- Dilute the solution to 10 mL with DI water, mix thoroughly, and transfer to clean and dry Erlenmeyer flask (sealed with a cork stopper)

	Sample No.	0.10 mg/mL $\text{Co}^{2+}$ (mL)	6 M HCl	50% KSCN	Acetone
Blank solution (step 2) →	0	0	0.8 mL	2.0 mL	4.8 mL
For calibration curve (step 3) {	1	0.50			
	2	1.00			
	3	1.50			
	4	2.00			
Unknown sample (step 4) →	Unknown	$0.5 \leq x \leq 2.0$			



# Step 5: Blank Adjust





Blank solution

- Use a dropping pipet to rinse the cuvette twice with blank solution, then add blank solution into the cuvette to  $\sim \frac{1}{2}$  height
- Hold only the top part of the cuvette, and clean off any residual liquid outside the cuvette with a lens tissue
- Place the cuvette gently into the cuvette holder and align the vertical white lines
- Close the lid, then press the **“BLANK”** button to blank adjust

✓ Do NOT use regular test tubes in place of cuvette



## Step 6: Absorbance Measurement



- Perform the measurement in the sequence of 1 → 4 (lower to higher concentration)
- Pour the solution back to the original flask, then rinse the cuvette with the next tested solution twice; filled to  $\frac{1}{2}$  height
- Wipe off any residual liquid outside the cuvette, then place the cuvette into the cuvette holder (align the vertical white lines)
- Close the lid, and read the absorbance value at 620 nm



# Clean-Up and Check-Out

- Wash all 50 mL Erlenmeyer flasks and place them in the oven
- Waste solution that contains  $\text{Co}^{2+}$  ions and organic solvent (acetone) should be disposed into the designated container in the fume hood
- Turn off the spectrophotometer and put the dust cover on
- Clean cuvette only with water (**to avoid scratching the cuvette, do not use brush to clean**), return the cuvette to TA
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- This is a **Full Report** experiment:
  - Member A: **Hand in prelab to the TA**
  - Member B: **Have the lab notes and results checked by the TA, and hand in the report next week**
- Groups on duty shall stay and help clean up the lab




# Plot the Calibration Curve

- 5-point scattered plot (including the blank solution) and label both the x and y axes properly; Include both the data table (Table 1) and the calibration curve (Figure 1) in your report
- Plot the linear regression line ( $y = mx + b$ ) – show the result and the  $R^2$  value together with the plot
- Calculate the original concentration of unknown  $\text{Co}^{2+}$  solution prior to dilution; the absorbance of the unknown sample should fall within the range of the calibration curve

**Table 1**

$[\text{Co}^{2+}]$ (mg/mL)	Abs.
0	0.000
0.0050	0.183
0.010	0.399
0.015	0.604
0.020	0.809

**Figure 1**





# Final Report (Full Version)

- Four experiments (E3, E8, E10, E12)
- Complete the data analysis and calculation part in the lab manual
- Plot data correctly and discuss potential sources of errors
- Hand in the report in the following week together with the prelab and lab records
- 50 points per report (5 pts deduction for late submission < 1 week)

## I. Prelab exercise

- ✓ Objectives
- ✓ Principles
- ✓ Chemicals
- ✓ Procedures

**15 pts**

**+**

## II. Lab Notes

- ✓ Observation
- ✓ Operation
- ✓ Reaction condition
- ✓ Data and results

**10 pts**

## III. Final report

- ✓ Data analysis
- ✓ Elaborate results
- ✓ Error analysis
- ~~✓ Questions and discussion~~

**+**

**25 pts**



# Lab Report Grading Rubrics

Category	Guidelines	Pts
I. Prelab exercise	1. Briefly summarize main principles and relevant equations	5
	2. List the chemicals' toxicity and physical and chemical properties	5
	3. Use flow chart to explain the experimental procedures	5
II. Lab notes	4. Record data with correct significant figures and units	5
	5. Record observations, operations, and reaction conditions in details	5
III. Final report	6. Process data correctly (calculation included)	5
	7. Present final results with correct significant figures and units	5
	<i>8. Analyze the results with appropriate error discussions*</i>	<i>5</i>
	<i>9. Plot the results with correct XY axes and labeling*</i>	<i>5</i>
	<i>10. Elaborate findings and provide constructive suggestions*</i>	<i>5</i>

*\*Only for full reports*






# T12.2 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 1

- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid to the 0.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Press valve **E** to drain the liquid to the 5.00 mL marking
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use






# T12.3 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 2


- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid until it rises above the 5.00 mL marking
- Remove the pipet filler and quickly use an index finger to close the top of pipet. Use the finger to adjust the liquid level to the 5.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Release the index finger so that liquid is transferred
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use





# T13 – Volumetric Flask


- Clean the volumetric flask thoroughly, then rinse it with a small amount of solvent
- Using a funnel, transfer the solution to be diluted into the volumetric flask
- Fill solvent into the flask until about half full, swirl the flask to let the solution mix
- Add more solvent so that the liquid level approaches (but does not exceed) the inscribed mark
- Use a dropper pipet to add solvent slowly, so that the liquid level matches the inscribed mark
- Install the stopper cap (hold with a finger), invert the flask several times to ensure thorough mixing
- Pour the solution into a beaker for later use (do not store solution in the flask)
- Wash the volumetric flask immediately after use and let it air dry (do not put flask on a hot plate or in an oven)





# T17 – Spectrophotometer

- Turn on the power switch and let the instrument warm up for at least 20 minutes
- Ensure the *cuvette holder* is empty
- Press the “**Mode**” button several times until “A” (absorbance) appears on the screen
- Set the wavelength to the desired value (e.g. 620 nm)
- Press the “**Blank**” button to zero the reading
- Place a *cuvette* with blank solution into the cuvette holder. Align the white line on the cuvette toward you (do NOT use regular test tubes in the spectrophotometer)
- Press the “**Blank**” button to calibrate
- Place a sample solution into cuvette holder
- Close the lid of the sample compartment, record the absorbance reading



*Note – the T17 video shows the older Spectronic 20 model instead of the currently used SP-830*



# Lab Dispenser

- 1) Check the pre-set volume on the dispenser. Do not change the setting unless instructed to do so
- 2) Place the receiving flask under the tip of dispenser
- 3) To remove the air bubbles in the dispenser, lightly pull the piston pump up and down for several times
- 4) Gently pull the piston pump up until it reaches the end of travel range, then slowly push the piston down to obtain the solution

