



General Chemistry Laboratory

The Molar Volume of Nitrogen Gas



Preparation

- Put on your lab coat and safety goggles
- Turn off your mobile phone
- Place your backpack in the drawer or the cabinet
- Put your prelab on lab bench (hold it down with something heavy) for ATA to sign

Collect the following items

- Two small test tubes (oven)
- An iron stand and two three prong clamps
- 250 mL Erlenmeyer flask, 500 mL Florence flask, a rubber stopper with rubber tube and glass tube (boxes on the central islands)




Objective and Principles

- **Objective:** Using the ideal gas law to deduce the molar volume of nitrogen gas at STP
- **Lab techniques:**
 - Using an analytical balance to weigh chemicals
 - Measuring volume, pressure, and temperature
- **Definitions:**
 - Molar volume: the volume of a mole of substance
 - STP (standard temperature and pressure): 0°C, '**1 atm**'
 - For ideal gas, molar volume = 22.414 L at STP

✓ STP: Abbreviation for standard temperature (273.15 K or 0°C) and pressure (105 Pa). Ref: *IUPAC. Compendium of Chemical Terminology, 2nd ed. (the "Gold Book")*



Nitrogen-Producing Reaction



Initial:	n_1	n_2			
Change:	$-n_1$	$-n_1$	$+n_1$	$+n_1$	$+n_1$
Final:	0	$n_2 - n_1$	n_1		n_1

- Sulfamic acid ($\text{H}_2\text{NSO}_3\text{H}$, n_1 mole) being the limiting reagent
- Sodium nitrite (NaNO_2) being the excess reagent ($n_2 > n_1$)



Experimental Setup

Fastened by a three prong clamp

Fastened by a three prong clamp
(Use the iron stand on hot plate)

Erlenmeyer flask
(A)

Small test tube

Florence flask
(B)


Rubber/glass
tube (C)

400 mL beaker
(D)

Lab bench



Nitrogen-Producing Reaction



- Sulfamic acid ($\text{H}_2\text{NSO}_3\text{H}$, n_1 mole) being the limiting reagent
- Sodium nitrite (NaNO_2) being the excess reagent ($n_2 > n_1$)
- As per ideal gas law, the molar volume of N_2 at STP (0°C , 1 atm) can be related to the volume of N_2 (V_1) at room temperature (T_1) given the number of mole (n_1):

$$\frac{1(\text{atm}) \times V_{STP}(\text{L})}{1(\text{mol}) \times 273.15(\text{K})} = \frac{P_{\text{N}_2} \times V_1}{n_1 \times T_1} = \frac{(P_{\text{atm}} - P_{\text{H}_2\text{O}}) \times \Delta V}{n_1 \times T_1}$$

Mercury barometer (points to P_{atm})
Vapor pressure of water (points to $P_{\text{H}_2\text{O}}$)
Volume of collected water (points to ΔV)
mole $\text{H}_2\text{NSO}_3\text{H}$ (points to n_1)
Thermometer (points to T_1)



Step 1/6: Measuring Chemicals

- Measure roughly 1 g NaNO_2 and dissolve them with 50 mL DI water in Erlenmeyer flask A
- Measure the weight of an empty small test tube (W_1) using an analytical balance
- Place roughly 1.0-1.1 g $\text{H}_2\text{NSO}_3\text{H}$ into the small test tube, record the accurate weight (W_2)



✓ Use a small beaker to hold the test tube → only the weight difference ($W_2 - W_1$) matters



Step 1/6: Measuring Chemicals



Place a small beaker
(container)



Close
windshield
& zeroing
(TARE)



Display shows 0.0000 g
(Deduction of container's mass)



Place test
tube in &
close
windshield



W_1 : 7.5757 g
(Mass of small test tube)

Take out
small test tube



Mass of $\text{H}_2\text{NSO}_3\text{H}$
 $W_2 - W_1 = 1.0445 \text{ g}$
(Five s.f.)



W_2 : 8.6202 g
(Mass of test tube and $\text{H}_2\text{NSO}_3\text{H}$)



Place test tube
in and close
windshield




Use the skinny end of a spatula to
put ~1 g of $\text{H}_2\text{NSO}_3\text{H}$ (ca. 5 times)
in the small test tube



Step 2/6: Check the Water Flow

- Fill both the Florence flask B and the beaker D with water
- Fill the rubber/glass tube C with water, apply the pinch clamp
- Fix the Florence flask with a three prong clamp to avoid breaking
- Use the rubber/glass tube C to connect B and D
- Loosen the pinch clamp and check if the water can flow freely between the two container; re-apply the clamp




- ✓ The end of glass tube shaft should nearly touch the bottom of the Florence flask
- ✓ Check whether there is air bubbles in the rubber/glass tube



Step 3/6: Set up the Reaction Vessel


- Use a three prong clamp to fix the Erlenmeyer flask A at a tilted angle
- Place the small test tube containing $\text{H}_2\text{NSO}_3\text{H}$ at the bottom of Erlenmeyer flask A using a tweezer
- Install the rubber stopper onto the Erlenmeyer flask

- ✓ Practice with the empty test tube first
- ✓ If NaNO_2 is accidentally mixed with $\text{H}_2\text{NSO}_3\text{H}$, clean the glassware and re-weight both chemicals





Step 4/6: Balance the Pressure




- Ensure all rubber stoppers are tightly fitted
- Loosen the pinch clamp, adjust the height of beaker so that the water levels in B & D become equal

- Tighten the pinch clamp and dispose all the remaining water in the beaker D
- Measure the weight of empty beaker





Ask a TA/ATA to Check Your Setup





Step 5/6: Start Generating Nitrogen

- Straighten the Erlenmeyer flask and let the chemicals mix
- Loosen the pinch clamp quickly
- Swirl the flask gently and observe water flow
- Maintain the rubber tube beneath the water level in the beaker

- ✓ Rubber tubes shouldn't be twisted
- ✓ Brown NO_2 gas may be produced via a side reaction





Step 6/6: Adjust Pressure and Temp.

- Place the Erlenmeyer flask A in a room-temperature water bath
- Adjust the height of the Florence flask so that the water levels in B & D become equal
- Re-apply the pinch clamp onto the rubber tube C
- Measure the weight of repelled water in beaker $\rightarrow \Delta V$
- Record room temperature T_1 and pressure P_{atm}
- Use Appendix 7 to find P_{H_2O} (vapor pressures of water)

$$\frac{1(atm) \times V_{STP}(L)}{1(mol) \times 273.15(K)} = \frac{(P_{atm} - P_{H_2O}) \times \Delta V}{n_1 \times T_1}$$

$$V_{STP} = ?$$



Clean-Up and Check-Out

- Clean and return the small test tubes
- Pour the solution waste into the drain directly
- Tuck the lab stools underneath the lab bench
- Clean up the lab bench and check personal equipment inventory (have an associate TA signed the check list)
- **Use the correct significant figures and units** (e.g. 1.0445 g, 25.13 °C, and 359.12 mL)
- This is a **Brief Report** experiment:
 - Complete calculation using correct significant figures
 - **Hand in prelab/lab note/report together to the TA**
- Groups on duty shall stay and help clean up the lab



Notes and Reminders

- Wear personal protective equipment (PPE, i.e. lab coat, safety goggles, closed-toe shoes, long pants) at all time in the laboratory
- Bring a scientific calculator (smartphone is not allowed)
- Communicate with your lab buddy
- Communicate with TA/ATA should you have any question



Electronic Balance

- Unless instructed, do not move the balance so that proper calibration is maintained
- Do not overload the balance (the maximum load is 610 grams for *electronic balance*, and 210 grams for *analytical balance*)
- Before use, warm up the balance for at least 30 min and ensure that it is level and clean
- Do not put chemicals directly on the weighing pan – use a folded weighing paper, a weighing boat, or a beaker (mind the weight limit) as container
- Close all windshields on the *analytical balance* before zeroing and recording values
- Maintain the tidiness inside and outside the balance; use the provided soft brush to clean spillages
- Do not weigh hot objects as the convective airflow will affect the measured mass



Electronic Balance (resolution 0.01 g)



Analytical Balance (resolution 0.0001 g)



Weighing Chemicals

- Read the label on the chemical bottle carefully before proceeding to weigh
- For solid chemicals, place a folded weighing paper or a beaker on the electronic balance to hold chemicals. Use a clean and dry spatula to move chemicals
- For liquid chemicals, use a clean and rinsed dropper pipet
- Unless specifically instructed, never return any excess chemical to the original bottle to avoid contamination – use the designated waste bin
- Maintain the tidiness inside and outside the balance – move appropriate amount with spatula to avoid any spillage, and use the provided soft brush to clean scattered chemicals
- Close the cap of chemical bottle immediately after use




Figure T10-1 Weighing chemicals



Mercury Barometer




Figure T3-1 Illustration of mercury barometer