



General Chemistry Laboratory

Quantitative Analysis of Vitamin C



Preparation

Collect the following items

- One 25 mL transfer pipet + one pipet filler
- Two 125 mL Erlenmeyer flasks
- One 25 mL buret + one glass funnel
- One 100 mL volumetric flask
- The TA will distribute one magnetic stir bar to each group

From your personal equipment

- Clean and dry one 100 mL beaker and place it in the oven (for taking KIO_3 solution later)
1. Accurately weigh two low-dosage (or one high-dosage) vitamin C tablets
 2. Wrap each tablet in weighing paper, then use a hammer to gently crush the tablet on a concrete surface (NOT on the lab bench)
 3. Add each portion of crushed tablet in 50 mL DI water (low-dosage: 125 mL Erlenmeyer flask; high-dosage: 100 mL beaker), stir to dissolve




Objective and Principles


- **Objective:** Determine the actual content of vitamin C in commercially available tablets and juices by reduction-oxidation titration
- **Lab techniques:**
 - Hot plate/magnetic stirrer
 - Lab dispenser
 - Volumetric flask
 - Transfer pipet
 - Performing titration using a buret



Redox Titration





- **Vitamin C** (L-Ascorbic acid, $C_6H_8O_6$) is a reducing agent, i.e. oxidized in the reaction
- I_2 (the oxidizing agent) is generated by the reaction between **KIO₃** and **NaI** in an acidic medium:



- A known concentration of KIO₃ (titrant) is used to titrate vitamin C solutions (concentration to be determined)
- Once the equivalence point (# mole vitamin C = # mole I₂) is passed, the additional I₂ would react with starch indicator to produce a **blue-black complex** that indicates the titration endpoint



Redox Titration



$$\frac{\text{IO}_3^-(\text{mol})}{1} = \frac{\text{I}_2(\text{mol})}{3} = \frac{\text{C}_6\text{H}_8\text{O}_6(\text{mol})}{3}$$

$$\text{C}_6\text{H}_8\text{O}_6 (\text{mol}) = \text{I}_2 (\text{mol}) = \text{IO}_3^- (\text{mol}) \times 3$$

Volume of KIO_3 solution used \rightarrow # mole $\text{IO}_3^- \rightarrow$ # mole vitamin C



Experiment Tasks

Vitamin C Tablet

Low-dosage (≤ 250 mg)

2 tablets

OR

High-dosage (> 250 mg)

1 tablet (2 titrations)

Juice

50 mL per portion

*** 2 portions**


Note: Avoid using foaming ingots or jelly drops (capsules are okay)

Note: Avoid carbonated drinks or juices of dark colors



Task 1-1: Low-Dosage Vitamin C Tablet

1. Accurately weigh two low-dosage vitamin C tablets
2. Wrap each tablet in weighing paper, then use a hammer to gently crush the tablet on a concrete surface (NOT on the lab bench)
3. Dissolve each crushed tablet with 50 mL DI water (graduated cylinder) in a 125 mL Erlenmeyer flask, then stir to dissolve



✓ Do NOT use the heating function of the hot plate



Task 1-1: Low-Dosage Vitamin C Tablet


4. For each vitamin C tablet sample solution, add 2 mL 1 M NaI, 2 mL 1 M HCl, and **1 mL 2% starch solution** (volumes are preset on the lab dispensers – do not adjust them)
5. Titrate the solution with **0.0250 M KIO₃** (read and record V_i and V_f to 0.01 mL) until the titration endpoint
6. Repeat the procedure on the other vitamin C tablet, and use the average result to calculate the ascorbic acid contents in each vitamin C tablet





Task 1-2: High-Dosage Vitamin C Tablet

1. Accurately weigh one high-dosage vitamin C tablet
2. Wrap the tablet in weighing paper, then use a hammer to gently crush the tablet on a concrete surface (NOT on the lab bench)
3. Dissolve the crushed tablet with 50 mL DI water in 100 mL beaker, then stir to dissolve
4. Transfer the solution to a 100 mL volumetric flask; rinse the beaker 2-3 times with DI water and combine the wash solution to the volumetric flask
5. Dilute the solution to 100 mL; transfer the diluted solution to another clean beaker



✓ Do NOT use the heating function of the hot plate



Task 1-2: High-Dosage Vitamin C Tablet

6. Use a 25 mL transfer pipet to take 25.0 mL of diluted sample solution to a 125 mL Erlenmeyer flask, add 2 mL 1 M NaI, 2 mL 1 M HCl, and **1 mL 2% starch solution** (volumes are preset on lab dispensers – do not adjust them)
7. Titrate the solution with **0.0250 M KIO₃** (read and record V_i and V_f to 0.01 mL) until the titration endpoint
8. Take another portion of 25.0 mL diluted sample solution and repeat the titration procedure; use the average result to calculate the ascorbic acid contents in the original vitamin C tablet





Task 2: Vitamin C Content in Juice

1. Use a 25 mL transfer pipet to take **50.0 mL** of package juice (do the transfer twice) into a 125 mL Erlenmeyer flask
2. Add 2 mL 1 M NaI, 2 mL 1 M HCl, and 1 mL 2% starch solution (volumes are preset on lab dispensers)
3. Titrate the solution with 0.0250 M KIO_3 (read and record V_i and V_f to 0.01 mL) until the titration end point
4. Take another portion of 50.0 mL juice sample and repeat the titration procedure; use the results to calculate the ascorbic acid contents in the package juice





Result Calculation



Task 1 – Vitamin C tablet (___ mg/tablet)

- Low-dosage: $(\Delta V \text{ mL} \times \underline{0.0250 \text{ M}}) \times \frac{3 \text{ C}_6\text{H}_8\text{O}_6}{1 \text{ IO}_3^-} \times 176.12 \frac{\text{mg}}{\text{mmol}}$
mmol IO₃⁻
- High-dosage: $(\Delta V \text{ mL} \times \underline{0.0250 \text{ M}}) \times \frac{3 \text{ C}_6\text{H}_8\text{O}_6}{1 \text{ IO}_3^-} \times 176.12 \frac{\text{mg}}{\text{mmol}} \times \frac{100.0 \text{ mL}}{25.0 \text{ mL}}$

Task 2 – Juice (___ mg/100 mL)

- $(\Delta V \text{ mL} \times \underline{0.0250 \text{ M}}) \times \frac{3 \text{ C}_6\text{H}_8\text{O}_6}{1 \text{ IO}_3^-} \times 176.12 \frac{\text{mg}}{\text{mmol}} \times \frac{100.0 \text{ mL}}{50.0 \text{ mL}}$




Additional Notes

- Use a clean and dry 100 mL beaker to take **0.0250 M standard KIO_3** solution (titrant) – any remaining water would affect the actual concentration of KIO_3
- Refer to slide for detailed instructions for using a transfer pipet
- Refer to slide for detailed instructions for using a volumetric flask
- Refer to slide for detailed instructions for titration
- Wash the volumetric flask thoroughly after use (do NOT put them into the oven)
- All iodine-containing waste solution should be collected and recycled
- Do not eat or drink during the lab session



Clean-Up and Check-Out


- Place the waste juice bottles into the designated recycle bin
- The cleaned buret should be clamped upside-down on the buret clamp
- 25 mL transfer pipet: rinse with water, then clean in a sonicator for 1 min, washed and clamped upside-down on the buret clamp
- Return the magnetic stir bar to TA
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- This is a **Brief Report** experiment:
 - Complete calculation using correct significant figures
 - **Hand in prelab/lab note/report together to the TA**
- Groups on duty shall stay and help clean up the lab





Transfer (Volumetric) Pipet


- Clean the pipet and rinse it twice with a small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of the pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid until it rises above the inscribed line
- Remove the pipet filler and quickly use an index finger to close the top of pipet
- Use a finger to adjust the liquid level to the inscribed line. Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Release the index finger so that the liquid is transferred
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use





Volumetric Flask




- Clean the volumetric flask thoroughly, then rinse it with a small amount of solvent
- Using a clean funnel, transfer the solution to be diluted into the volumetric flask
- Fill solvent into the flask until it is about half full, and swirl the flask to let the solution mix
- Add more solvent so that the liquid level approaches (but does not exceed) the inscribed mark
- Use a dropper pipet to add solvent slowly so that the liquid level matches the inscribed mark
- Install the stopper cap (hold with a finger), invert the flask several times to ensure thorough mixing
- Pour the solution into a beaker for later use (do not store the solution in the flask)
- Wash the volumetric flask immediately after use and let it air dry (do not put the flask on a hot plate or in an oven)





Titration

- Clean the buret with DI water, then rinse twice with ~5 mL of titrant (use a funnel to add titrant)
- Open the stopcock to repel the air at the buret tip
- Adjust the height of the buret so that its tip is slightly lower than the lid of the receiving flask
- Read and record the initial volume (V_i) on the buret to 0.01 mL
- With the stopcock on the right side, use your non-dominant hand to control the stopcock while the dominant hand swirls the receiving flask in a circular motion
- At the titration endpoint, read and record the final volume (V_f) on the buret to 0.01 mL
- After the experiment, wash the buret and let it dry upside-down on the buret clamp





Lab Dispenser

- 1) Check the pre-set volume on the dispenser. Do not change the setting unless instructed to do so
- 2) Place the receiving flask under the tip of the dispenser
- 3) To remove the air bubbles in the dispenser, lightly pull the piston pump up and down several times
- 4) Gently pull the piston pump up until it reaches the end of the travel range, then slowly push the piston down to obtain the solution

