



General Chemistry Laboratory

Iodine Clock – The Initial Rate Method



Preparation

Collect the following items

- Ten 50 mL Erlenmeyer flasks (clean and oven dry)
- Two cork stoppers
- Two 5 mL graduated pipets and one pipet filler
- One stop watch (distributed by TA)

From your personal equipment

- Two 100 mL beakers (clean and oven dry) – labelled with $K_2S_2O_8$ and K_2SO_4 , respectively
- A scientific calculator




Objective and Principles

- **Objective:** Determine the *rate law* of the chemical reaction between persulfate ($\text{S}_2\text{O}_8^{2-}$) and iodide (I^-) ions




- A limiting amount of thiosulfate ion ($\text{S}_2\text{O}_3^{2-}$) is added as a measuring tool in order to determine the rate of the above reaction:




- **Lab techniques:**
 - Graduated pipet
 - Lab dispenser
 - Measuring initial reaction rate



Determine the Rate of a Reaction



Rate and rate law to be determined



$\text{S}_2\text{O}_3^{2-}$ being the limiting reagent


- The second reaction occurs much faster than the first reaction
→ I_2 formed by the first reaction is consumed immediately by the second reaction
- When the limiting reagent $\text{S}_2\text{O}_3^{2-}$ is used up, I_2 starts to accumulate and combine with I^- to form I_3^- , which shows purple-blue color in the presence of starch indicator

$$\Delta[\text{S}_2\text{O}_3^{2-}] = 2 \times \Delta[\text{S}_2\text{O}_8^{2-}]$$

$$\text{rate} = \frac{-\Delta[\text{S}_2\text{O}_8^{2-}]}{\Delta t} = \frac{-\frac{1}{2}\Delta[\text{S}_2\text{O}_3^{2-}]}{\Delta t}$$



Initial Rate Method



$$\text{Rate} = k[\text{S}_2\text{O}_8^{2-}]^m[\text{I}^-]^n$$

Table 1 Volumes of reagents for the initial rate method (total volume of 10.0 mL)

Trial No.	0.20 M NaI (mL)	0.20 M NaCl* (mL)	0.0050 M Na ₂ S ₂ O ₃ (mL)	2% Starch (mL)	0.10 M K ₂ SO ₄ * (mL)	0.10 M K₂S₂O₈ (mL)	Reaction time Δt (s)
1	2.0	2.0	1.0	1.0	2.0	2.0	109
2	2.0	2.0	1.0	1.0	0	4.0	59
3	4.0	0	1.0	1.0	2.0	2.0	58

$[\text{I}^-] \uparrow$
 $\Delta t \downarrow$
rate \uparrow


Limiting reagent
Fixed # of moles

$[\text{S}_2\text{O}_8^{2-}] \uparrow$
 $\Delta t \downarrow$
rate \uparrow

*NaCl and K₂SO₄ are added to maintain the ionic strength of the solution



Reaction Order



$$\text{Rate} = k[\text{S}_2\text{O}_8^{2-}]^m[\text{I}^-]^n$$

$$r = \frac{-\Delta[\text{S}_2\text{O}_8^{2-}]}{\Delta t} = \frac{-\frac{1}{2}\Delta[\text{S}_2\text{O}_3^{2-}]}{\Delta t} = \frac{-\frac{1}{2}(0 - 0.00050)}{\Delta t} = \frac{0.00025}{\Delta t}$$

Rate is inversely proportional to Δt

$$\frac{\text{rate}_2}{\text{rate}_1} = \frac{0.00025/\Delta t_2}{0.00025/\Delta t_1} = \frac{\Delta t_1}{\Delta t_2} = \frac{109}{59} = \frac{k(2.0[\text{S}_2\text{O}_8^{2-}]_1)^m([\text{I}^-]_1)^n}{k([\text{S}_2\text{O}_8^{2-}]_1)^m([\text{I}^-]_1)^n} = (2.0)^m$$

$$\frac{\text{rate}_3}{\text{rate}_1} = \frac{0.00025/\Delta t_3}{0.00025/\Delta t_1} = \frac{\Delta t_1}{\Delta t_3} = \frac{109}{58} = \frac{k([\text{S}_2\text{O}_8^{2-}]_1)^m(2.0[\text{I}^-]_1)^n}{k([\text{S}_2\text{O}_8^{2-}]_1)^m([\text{I}^-]_1)^n} = (2.0)^n$$

The rate constant k can be calculated with m and n known



Step 1: Prepare Sample Solutions

- Wash ten 50 mL Erlenmeyer flasks, dry them in oven and **let cool**
- Use 100 mL beakers to take **~30 mL** 0.10 M $K_2S_2O_8$ and 0.10 M K_2SO_4
- Reference to Table 1, use lab dispensers to measure 0.20 M NaI, 0.20 M NaCl, 0.0050 M $Na_2S_2O_3$, and 2% starch solutions into the first two Erlenmeyer flasks (two copies for each entry)
- Use a 5 mL graduated pipet to measure 0.10 M K_2SO_4 accurately, then transfer to Erlenmeyer flasks




- ✓ Use the same set of chemicals through out this lab
- ✓ Rinse graduate pipets twice before use



Step 2: Measure the Reaction Time

- Use a 5 mL graduated pipet to measure the required amount of 0.10 M $K_2S_2O_8$, add into the Erlenmeyer flask and start the timer immediately
- Install the cork stopper and swirl the Erlenmeyer flask for 20 s
- Leave the Erlenmeyer flask on the benchtop and record the time it takes for the purple-blue color to appear




✓ Remove the pipet filler to expel the $K_2S_2O_8$ solution from the graduated pipet



Step 3: Calculate Reaction Orders

- Repeat each trial twice (if the reaction times Δt differ for more than 3 s, redo the trial one more time)
- Calculate the average Δt values
- Calculate the values of m and n , then calculate the rate constant k



$$\text{Rate} = k[\text{S}_2\text{O}_8^{2-}]^m[\text{I}^-]^n$$

✓ Keep two significant figures for reaction orders m and n

Example:

Reaction time Δt (s)			
Trial 1	1'48''	1'54''	1'49''
Trail 2	58	59	
Trail 3	53	57	59

$$\text{Average } \Delta t_1 = 108.5 = 109 \text{ (s)}$$


$$\text{Average } \Delta t_2 = 58.5 = 59 \text{ (s)}$$

$$\text{Average } \Delta t_3 = 58 \text{ (s)}$$



Step 4: Participate in Iodine-Clock Symphony Competition

- Obtain an assigned reaction time Δt from the TA, and design a set of reaction conditions to achieve it (the assigned times will allow the color change of solution to match the beats in the 'Habanera' aria of opera 'Carmen')
- Verify the reaction conditions on your own, and make adjustment as needed
- At the scheduled competition time, bring your Erlenmeyer flask (containing all reactants except for $\text{K}_2\text{S}_2\text{O}_8$) and graduated pipet (containing $\text{K}_2\text{S}_2\text{O}_8$) to the podium
- On TA's mark, add $\text{K}_2\text{S}_2\text{O}_8$ into the Erlenmeyer flask, swirl for 20 s, then place the Erlenmeyer flask in the designated area



- ✓ The total volume of solution should be fixed at 10 mL
- ✓ Balance the ionic strength with NaCl and K_2SO_4



How to Design Reaction Conditions

Example 1: Target time is 90 s (if m has been determined to be **0.90**)

Trial No.	0.20 M NaI (mL)	0.20 M NaCl (mL)	0.0050 M Na ₂ S ₂ O ₃ (mL)	2% Starch (mL)	0.10 M K ₂ SO ₄ (mL)	0.10 M K₂S₂O₈ (mL)	Reaction time Δt (s)
1	2.0	2.0	1.0	1.0	2.0	2.0	109
Rate ₉₀	2.0	2.0	1.0	1.0	4.0 - x	x	90

By dispenser
(do not adjust the pre-set volume !)

By graduated pipet

$$\frac{\text{rate}_1}{\text{rate}_{90}} = \frac{\Delta t_{90}}{\Delta t_1} = \frac{90}{109} = \frac{k[\text{S}_2\text{O}_8^{2-}]^{0.90} [\text{I}^-]^{0.90}}{k[\text{S}_2\text{O}_8^{2-}]^{0.90} [\text{I}^-]^{0.90}} = \left(\frac{2.0}{x}\right)^{0.90}$$

$$\log\left(\frac{90}{109}\right) = 0.90 \log\left(\frac{2.0}{x}\right)$$

$$x = 2.47 = 2.5 \text{ mL } (x < 4.0 \text{ mL})$$



How to Design Reaction Conditions

Example 2: Target time is 40 s (if m has been determined to be **0.90**)

Trial No.	0.20 M NaI (mL)	0.20 M NaCl (mL)	0.0050 M Na ₂ S ₂ O ₃ (mL)	2% Starch (mL)	0.10 M K ₂ SO ₄ (mL)	0.10 M K₂S₂O₈ (mL)	Reaction time Δt (s)
3	4.0	0	1.0	1.0	2.0	2.0	58
Rate ₄₀	4.0	0	1.0	1.0	4.0 - y	y	40

By dispenser (do not adjust the pre-set volume !)

By graduated pipet

$$\frac{\text{rate}_3}{\text{rate}_{40}} = \frac{\Delta t_{40}}{\Delta t_3} = \frac{40}{58} = \frac{k[\text{S}_2\text{O}_8^{2-}]^{0.90}[\text{I}^-]^{0.90}}{k[\text{S}_2\text{O}_8^{2-}]'^{0.90}[\text{I}^-]^{0.90}} = \left(\frac{2.0}{y}\right)^{0.90}$$

$$\log(40/58) = 0.90 \log(2.0/y) \quad \boxed{y = 3.02 = 3.0 \text{ mL } (y < 4.0 \text{ mL})}$$



Clean-Up and Check-Out


- Wash all 50 mL Erlenmeyer flasks and place them in the oven
- Waste solution containing iodine should be disposed into the designated container
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- This is a **Full Report** experiment:
 - Member A: **Have the lab notes and results checked by the TA, and hand in the report next week**
 - Member B: **Hand in prelab to the TA**
- Groups on duty shall stay and help clean up the lab



T12.2 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 1

- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid to the 0.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Press valve **E** to drain the liquid to the 5.00 mL marking
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use






T12.3 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 2


- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid until it rises above the 5.00 mL marking
- Remove the pipet filler and quickly use an index finger to close the top of pipet. Use the finger to adjust the liquid level to the 5.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Release the index finger so that liquid is transferred
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use





Lab Dispenser

- 1) Check the pre-set volume on the dispenser. Do not change the setting unless instructed to do so
- 2) Place the receiving flask under the tip of dispenser
- 3) To remove the air bubbles in the dispenser, lightly pull the piston pump up and down for several times
- 4) Gently pull the piston pump up until it reaches the end of travel range, then slowly push the piston down to obtain the solution





Final Report (Full Version)

- Four experiments (E5, E8, E10, E12)
- Complete the data analysis and calculation part in the lab manual
- Plot data correctly and discuss potential sources of errors
- Hand in the report in the following week together with the prelab and lab records
- 50 points per report (5 pts deduction for late submission < 1 week)

I. Prelab exercise

- ✓ Objectives
- ✓ Principles
- ✓ Chemicals
- ✓ Procedures

II. Lab Notes

- ✓ Observation
- ✓ Operation
- ✓ Reaction condition
- ✓ Data and results

III. Final report

- ✓ Data analysis
- ✓ Elaborate results
- ✓ Conclusion
- ✓ Error analysis

15 pts

+

10 pts

+

25 pts