



# General Chemistry Laboratory

---

## Buffer Solutions



# Preparation

## Collect the following items

- One 250 mL volumetric flask
- Two 10 mL graduated pipets and one pipet filler
- Four 50 mL beakers
- Two plastic drop pipets**
- The TA will distribute one stir bar to each group

## From your personal equipment

- Two test tubes (clean and dry), one 400 mL beaker
- 10 mL and 50 mL graduated cylinders
- One wash bottle and one 1 L plastic beaker



# Objective and Principles

## ■ Objective:

- Prepare HOAc/NaOAc buffer solutions
- Investigate the factors affecting buffer capacity: **concentration** and  **$[A^-]/[HA]$  ratio**


## ■ Lab techniques:

- Operate a pH meter
- Volumetric flask
- Graduated pipet



# Buffer Solution

- Function: Resisting pH change upon the addition of either acid or base
- Composition:
  - A mixture of weak acid + its conjugate base (e.g. HF/NaF)
  - A mixture of weak base + its conjugate acid (e.g. NH<sub>3</sub>/NH<sub>4</sub>Cl)
- Acid-base equilibrium:



$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

- Henderson-Hasselbalch equation:

$$[\text{H}^+] = \frac{K_a [\text{HA}]}{[\text{A}^-]} \longrightarrow \text{pH} = -\log[\text{H}^+] = -\log K_a - \log\left(\frac{[\text{HA}]}{[\text{A}^-]}\right)$$
$$\text{pH} = \text{p}K_a + \log\frac{[\text{A}^-]}{[\text{HA}]}$$



# Buffer Capacity

- Definition: The amount of  $\text{OH}^-$  or  $\text{H}^+$  (in moles) required to change the pH of buffer solution by 1 unit
- A buffer solution would have a higher buffer capacity when:
  - The total concentration of two components ( $[\text{HA}] + [\text{A}^-]$ ) is higher
  - The concentration ratio  $[\text{A}^-]/[\text{HA}]$  is closer to 1
- **The most effective buffer range:  $\text{p}K_a \pm 1$** 
  - When  $[\text{A}^-] = 10 \times [\text{HA}]$ ,  $\text{pH} = \text{p}K_a + 1$
  - When  $[\text{A}^-] = [\text{HA}]/10$ ,  $\text{pH} = \text{p}K_a - 1$



# Working Principles of pH Meter

A pH meter consists of three parts:

➤ **pH electrode assembly**

- Reference electrode (often Ag/AgCl) – potential is fixed
- Indicator electrode (glass frit) – potential varies with  $[H^+]$


➤ **Voltmeter**: measure the potential difference ( $E_m$ ) between the two electrodes

➤ **Thermoprobe**: measure the temperature of solution





# Working Principles of pH Meter



$$E_m = \underbrace{K}_{\text{Reference electrode}} - \underbrace{2.3RT(\text{pH})/nF}_{\text{Indicator electrode}}$$

*Reference  
electrode*

*Indicator  
electrode*

- $E_m$ : measured cell potential
- $K$ : constant, determined by the type of electrode used
- $R$ : gas constant
- $T$ : absolute temperature of the solution
- pH: pH value of solution
- $n$ : number of moles of electrons transferred in the reaction
- $F$ : Faraday constant

- The pH meter needs to be calibrated in standard solutions (pH 7.00 and pH 4.00)
- After calibration, the measured  $E_m$  can be converted to pH

$$E_m = K + mT \times (\text{pH})$$



# Experiment Tasks


1. Calibrate pH-meter



2. Prepare HOAc and NaOAc solution



3. Measure the volume of one drop of HCl/NaOH



4. Prepare buffer solutions and record their pH values



50 mL beaker

Add 5 drops of HCl(aq)

Add 5 drops of NaOH(aq)

Record the pH values dropwise





# Step 1: Calibration of pH Meter



- Press the **“POWER”** button to turn on the pH meter. Let it warm up for at least 10 min
- Press the **“HOLD”** button to suspend pH reading
- Remove the electrode cap from by rotating it
- Use a wash bottle to rinse the electrode assembly, then wipe it dry gently with a tissue




✓ Each pH electrode assembly costs NTD 3,000 → Be careful when operate it



# Step 1: Calibration of pH Meter

- Press the “**MODE**” button several times until “Temp” appears on the screen. Check whether the temperature reading is close to RT
- Press the “**MODE**” button again to switch to “pH” function



- Immerse both thermoprobe and pH electrode into **pH 7.00** buffer solution
- Adjust **Calib** knob until ‘7.00’ is shown
- Rinse the thermoprobe and pH electrode with DI water
- Switch to **pH 4.00** buffer solution
- Adjust **Slope** knob until ‘4.00’ is shown



## Step 2: Prepare HOAc & NaOAc Solutions

### Odd Group #'s

#### Prepare 0.050 M HOAc

- Measure **0.74 mL** of conc. (17 M) HOAc
- Transfer to a 250 mL volumetric flask that contains some DI water
- Dilute to 250 mL
- pH of solution should be ~3



### Even Group #'s

#### Prepare 0.050 M NaOAc


- Measure **1.70 g** of NaOAc·3H<sub>2</sub>O, dissolve them in ~ 100 mL DI water in a beaker
- Transfer to a 250 mL volumetric flask (rinse the beaker 2-3 times)
- Dilute to 250 mL
- pH of solution should be ~8

- ✓ Two groups share one set of HOAc & NaOAc solutions
- ✓ Transfer the prepared solution into 400 mL beakers (labelled properly)



# Step 3: Measure Drop Volumes

- Measure 5 mL of 1.0 M HCl and 1.0 M NaOH separately into two dry test tubes
- Use a plastic drop pipet to add 30~50 drops of 1.0 M HCl into a 10 mL graduated cylinder; record the volume of solution
- Use another clean plastic drop pipet to add 30~50 drops of 1.0 M NaOH into a 10 mL graduated cylinder; record the volume of solution




1.0 M NaOH

1.0 M HCl

$$V_{\text{avg, HCl}} = \text{_____ mL/drop}$$

$$V_{\text{avg, NaOH}} = \text{_____ mL/drop}$$





# Step 4: Preparing Test Solutions

- Use 10 mL graduate pipets to measure volume accurately for HOAc and NaOAc solutions
- Make two copies of each test solution in 50 mL beakers (for adding acid and base, respectively)
- Perform the measurements in the sequence of **(a)**, **(b)**, (c), (d), then **(e)**


Test solution	0.050 M HOAc (mL)	0.050 M NaOAc (mL)	DI water (mL)
<b>(a)</b>	<b>15</b>	<b>15</b>	<b>0</b>
<b>(b)</b>	<b>5</b>	<b>5</b>	<b>20</b>
(c)	30	0	0
(d)	0	30	0
<b>(e)</b>	<b>0</b>	<b>0</b>	<b>30</b>



# Step 5: Measuring pH Changes



Stir bar  
50 mL beaker



Salt  
bridge

- Both the thermoprobe and the pH electrode should be fixed on the holder and immersed in the test solution (the salt bridge of pH electrode should be fully immersed)
- Adjust the position of beaker and the pH electrode so that the stir bar does not collide with the electrode assembly
- Use plastic drop pipets to add HCl or NaOH into the test solution, record the solution pH dropwise
- When changing between solutions: press the “HOLD” button on the pH meter, rinse both the electrode and the thermoprobe with DI water, then wipe dry with a tissue paper



# Step 5: Measuring pH Changes



Stir bar  
50 mL beaker

Obtain **1.0 M HCl** with a drop pipet



Add one drop into soln. **(a)**



Mix thoroughly  
and record pH values



Repeat adding acid for 5 times



Prepare another part of soln. **(a)**



Obtain **1.0 M NaOH** with a drop pipet



Record pH values after adding each  
drop (5 drops in total)



Repeat the same procedures with soln.  
**(b)(c)(d)(e)**



# Clean-Up and Check-Out

- Immerse the pH electrode in DI water (use a 100 mL beaker)
- Switch off the pH meter (keep the power cord plugged in)
- Salt solutions resulted from acid-base neutralization can be disposed into the sink
- Return the magnetic stir bar to TA
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- This is a **Full Report** experiment:
  - Member A: **Have the lab notes and results checked by the TA, and hand in the report next week.**
  - Member B: **Hand in prelab to the TA**
- Groups on duty shall stay and help clean up the lab





# Example of Data Table


<b>Table 1</b>		<b>c</b>	<b>d</b>	<b>c</b>	<b>d</b>	<b>e</b>
	No. of drop	0.050 M HOAc/NaOAc	0.017 M HOAc/NaOAc	0.050 M HOAc	0.050 M NaOAc	DI water
Add NaOH(aq)	5	4.80	5.09	3.65	11.51	11.49
	4	4.76	4.99	3.59	11.39	11.40
	3	4.73	4.89	3.48	11.22	11.27
	2	4.70	4.80	3.33	10.93	11.09
	1	4.67	4.73	3.15	10.78	10.79
Avg. pH of acid and base	<b>0</b>	<b>4.65</b>	<b>4.64</b>	<b>3.00</b>	<b>7.71</b>	<b>5.49</b>
Add HCl(aq)	-1	4.61	4.56	2.82	6.37	3.01
	-2	4.58	4.46	2.71	6.07	2.71
	-3	4.55	4.36	2.59	5.85	2.55
	-4	4.51	4.26	2.47	5.72	2.43
	-5	4.47	4.15	2.39	5.60	2.35



# Plot pH vs. # Drops of Acid/Base

- Use 'x-y scattering plot' function in MS Excel (or similar software); use lines to connect adjacent data points
- Plot all five dataset in a same figure (with appropriate axis labels and legends)
- The data table and figure should be included in the full lab report

Figure 1






# Calculating Buffer Capacity

- For solutions (a) and (b), use the central 5 points in the dataset (see example Table 2)
- Use *linear regression analysis* to obtain the **slope** and  **$R^2$  value**
- The data table and figure should be included in the full lab report

**Table 2**

No.	a	b
Drops	0.050 M HOAc/NaOAc	0.017 M HOAc/NaOAc
-2	4.58	4.46
-1	4.61	4.56
0	<b>4.65</b>	<b>4.64</b>
1	4.67	4.73
2	4.70	4.80

**Figure 2**





# Calculating Buffer Capacity

Buffer Capacity (meq/L · pH)

$$= \frac{1}{\text{slope}} (\text{drop/pH}) \times V_{\text{drop}} (\text{mL/drop}) \times C_M (\text{eq/L}) \times \frac{1000 (\text{mL/L})}{30 \text{ mL}}$$

$(V_{\text{avg. of HCl}} + V_{\text{avg. of NaOH}}) / 2$


0.050 M (1:1): 22 meq/L·pH

0.017 M (1:1): 7.8 meq/L·pH

**Table 2**

No.	a	b
Drops	0.050 M HOAc/NaOAc	0.017 M HOAc/NaOAc
-2	4.58	4.46
-1	4.61	4.56
0	<b>4.65</b>	<b>4.64</b>
1	4.67	4.73
2	4.70	4.80

**Figure 2**






# T12.2 – Measuring (Graduated) Pipet

Deliver 5 mL solution – Method 1

- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid to the 0.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Press valve **E** to drain the liquid to the 5.00 mL marking
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use






# T12.3 – Measuring (Graduated) Pipet

## Deliver 5 mL solution – Method 2


- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid until it rises above the 5.00 mL marking
- Remove the pipet filler and quickly use an index finger to close the top of pipet. Use the finger to adjust the liquid level to the 5.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Release the index finger so that liquid is transferred
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use





# T13 – Volumetric Flask

- Clean the volumetric flask thoroughly, then rinse it with a small amount of solvent
- Using a funnel, transfer the solution to be diluted into the volumetric flask
- Fill solvent into the flask until about half full, swirl the flask to let the solution mix
- Add more solvent so that the liquid level approaches (but does not exceed) the inscribed mark
- Use a dropper pipet to add solvent slowly, so that the liquid level matches the inscribed mark
- Install the stopper cap (hold with a finger), invert the flask several times to ensure thorough mixing
- Pour the solution into a beaker for later use (do not store solution in the flask)
- Wash the volumetric flask immediately after use and let it air dry (do not put flask on a hot plate or in an oven)






# T16.1 – pH Meter

## Preparation and calibration:

- Press the **“POWER”** button to turn on the pH meter. Let it warm up for at least 10 min
- Press the **“MODE”** button several times until **“Temp”** appears on the screen. Check whether the temperature reading is close to RT
- Press the **“MODE”** button again to switch to **“pH”** function
- Press the **“HOLD”** button to suspend pH reading
- Rotate the electrode cap (3 M KCl) to remove it
- Use a wash bottle with DI water to rinse the electrode, then gently dry it with a tissue paper
- Immerse both the electrode and the thermoprobe into the pH 7.00 buffer solution
- Press **“HOLD”** again to unfreeze the reading. Adjust the **“Calib”** knob until 7.00 is shown
- Press **“HOLD”** to suspend pH reading. Remove the electrode and the thermoprobe, rinse them with DI water then immerse them into the pH 4.00 buffer solution
- Press **“HOLD”** again to unfreeze the reading. Adjust the **“Slope”** knob until 4.00 is shown









# T16.2 – pH Meter

## Measuring pH:

- Use a beaker to hold the solution to be tested. Add a magnetic stir bar and place the beaker atop a stirrer
- When changing between solutions: press the “HOLD” button on the pH meter, rinse both the electrode and the thermoprobe with DI water, then wipe dry with a tissue paper
- The *salt bridge* at the bottom of electrode should stay below the liquid level, however avoid hitting the electrode with stir bar
- Immerse the electrode in clean DI water when not in use
- For long-term storage, put the electrode back into the 3 M KCl solution in the electrode cap



**Salt  
bridge**