



General Chemistry Laboratory

Potentiometric Titration of Acid-Base



Preparation

Collect the following items

- One magnetic stir bar (TA distribute)
- One glass funnel
- One 25 mL buret
- One 100 mL volumetric flask
- Two 125 mL Erlenmeyer flasks (check if broken)
- One 5 mL pipet and pipet filler (shared)
- pH 7.00 and pH 4.00 standard buffer solution (shared)

From your personal equipment

- Two 100 mL beakers
- One wash bottle and one 1 L plastic beaker

Note:

Measure 2.5 mL vinegar to titrate



Objective and Principles

■ Objective:

- To **prepare** and to **standardize** secondary-standard solutions
- To determine the **equivalence point** and concentration of acetic acid by using the **electric potential method**
- To determine the dissociation constant of acetic acid, K_a



■ Lab techniques:

- Operate a **pH meter**
- Use of **volumetric flask**, and **graduated pipet**
- Determine the equivalence point by titration curves



Standardization of Acid or Base

- **Primary standard:** substance with high purity and high molar mass
 - Common base: sodium carbonate (Na_2CO_3)
 - Common acid: potassium hydrogen phthalate (**KHP**)
- **Secondary standard:** standardized acid or base
 - Common base: sodium hydroxide (NaOH)
 - Common acid: hydrogen chloride (HCl)
- KHP is a **monoprotic weak** acid with structure
- The neutralization reaction of KHP with NaOH in a **1:1** stoichiometric ratio




$$C_{\text{NaOH}} \cdot V_{\text{NaOH}} = n_{\text{KHP}} = \frac{\text{Mass}_{\text{KHP}}}{204.22}$$



Acid-base Indicator

- Acid-base indicator: a weak organic acid or base
- Weak acid (HIn) and its conjugate base (In⁻) with different colors



$$K_a = \frac{[\text{H}_3\text{O}^+][\text{In}^-]}{[\text{HIn}]}$$


$$\text{pH} = \text{pK}_a + \log \frac{[\text{In}^-]}{[\text{HIn}]}$$

- According to the pH range of the equivalence point, choose the appropriate indicator to match the **end-point** with the **equivalence point**
- At the equivalence point, pH of the solution:
 - Strong acid/weak base titration: pH < 7
 - Weak acid/strong base titration: pH > 7
 - Strong acid/strong base titration: pH = 7



Equivalence Point

- The **pH** of the reacting solution **changes significantly** near the **equivalence point**
- Base on the color change of the acid-base indicator or monitoring the change in pH values to determine the **end-point**



Indicator	Acid form	pH range	Basic form
Methyl orange	Red	4.2~6.0	Orange
Bromothymol blue	Yellow	6.0~7.6	Blue
Phenolphthalein	Colorless	8.3~10	Pink red

Weak acid / strong base titration curve



Equivalence Point

1. Acid-base titration curve




The point on the curve with the maximum slope is the equivalence point

2. First derivative of titration curve

The maximum point is the equivalence point

3. Second derivative of the titration curve


X-intercept of line A-B is the equivalence point






Dissociation Constant of Weak Acid

- Acid-base neutralization reaction:



- Dissociation of weak acid




$$K_a = \frac{[\text{A}^-][\text{H}_3\text{O}^+]}{[\text{HA}]}$$

- At half-equivalence point

$$[\text{HA}] = [\text{A}^-]; [\text{H}_3\text{O}^+] = K_a$$

$$\text{pH} = \text{p}K_a$$

Weak acid – strong base titration curve





Dissociation Constant of Weak Acid

For example:

Equivalence volume = 37.50 mL

- Half-equivalence volume = 18.75 mL

$V = 18.00 \text{ mL}, \text{pH} = 4.60$


$V = 19.10 \text{ mL}, \text{pH} = 4.65$

- pH of the half-equivalence point = 4.63

$$K_a = \frac{[A^-][H_3O^+]}{[HA]}$$

- $\text{p}K_a = \text{pH} = 4.63, K_a = 2.3 \times 10^{-5}$

Weak acid – strong base
titration curve





Working Principles of pH Meter

A pH meter consists of three parts:

➤ **pH electrode assembly**

- Reference electrode (often Ag/AgCl) – potential is fixed
- Indicator electrode (glass frit) – potential varies with $[H^+]$

➤ **Voltmeter**: measure the potential difference (E_m) between the two electrodes

➤ **Thermoprobe**: measure the temperature of solution





Fig. T16-1 pH meter in lab (SUNTEX SP-701)





Working Principles of pH Meter



$$E_m = \underbrace{K}_{\text{Reference electrode}} - \underbrace{2.3RT(pH)/nF}_{\text{Indicator electrode}}$$

Reference
electrode

Indicator
electrode

$$E_m = K + mT \times (pH)$$

- The pH meter needs to be **calibrated** in standard solutions (pH 7.00 and pH 4.00)
- After calibration, the measured E_m can be converted to pH

- E_m : measured cell potential
- K : constant, determined by the type of electrode used
- R : gas constant
- T : absolute temperature of the solution
- pH: pH value of solution
- n : number of moles of electrons transferred in the reaction
- F : Faraday constant



Experiment Tasks



1. Prepare NaOH(aq)
2. Standardize NaOH(aq)
3. Calibrate pH-meter
4. Titrate vinegar with NaOH(aq)



Step 1: Prepare 0.1 M NaOH Solution




- Measure 10 mL of 1 M NaOH to a 100 mL volumetric flask that **contains some DI water**
- Dilute to 100 mL

- Stopper and invert the volumetric flask several times to mix thoroughly
- Transfer to a 100 mL beaker (label properly)

- Rinse buret twice with **5 mL** NaOH(aq)
- Fill the buret
- Read initial volume of buret (V_i) to 0.01 mL



Step 2: Standardize NaOH Solution



1. Dissolve with 50 mL
DI water

2. Add 2 d of phenolphthalein
3. Titrate with 0.1 M NaOH




- Measure ca. **0.2~0.22 g** KHP with analytical balance
- Record accurate weight
- Place into a 125 mL Erlenmeyer flask

- Titrate the solution to appear pink and persist for 30 s
- Record V_i and V_f (to **0.01 mL**)
- Carry out a duplicate test
- Calculate average concentration of NaOH

✓ KHP dissolves slowly



Step 3: Calibration of pH Meter



- Press the **“POWER”** button to turn on the pH meter. Let it **warm up** for at least **10 min**
- Press the **“HOLD”** button to suspend pH reading
- Remove the electrode cap by rotating it
- Use a wash bottle to rinse the electrode assembly, then wipe it dry gently with a tissue





- ✓ Remove electrode cap and operate with care
- ✓ Each pH electrode assembly costs NTD 3000



Step 3: Calibration of pH Meter




- Press the **“MODE”** button several times until **“Temp”** appears on the screen. Check whether the temperature reading is close to RT
- Press the **“MODE”** button again to switch to **“pH”** function



- Immerse both thermoprobe and pH electrode into **pH 7.00** buffer solution
- Adjust **Calib** knob until **‘7.00’** is shown
- Rinse the thermoprobe and pH electrode with DI water
- Switch to **pH 4.00** buffer solution
- Adjust **Slope** knob until **‘4.00’** is shown



Step 4: Titration of Vinegar



- Transfer **2.5 mL** vinegar into a 100 mL beaker
- Record the **brand** and **concentration** of vinegar
- Add 40 mL of DI water
- Add 2 d. of phenolphthalein
- Place stirring bar, the electrode, and thermoprobe in solution
- Titrate with standardized 0.1 M NaOH and **record V_i , V_f , and pH value** after each addition
 - At pH **< 5.5**: add ~ 1 mL of NaOH
 - At pH **5.5~10**: add ~ **0.2 mL** NaOH
 - At pH **10~11**: add ~ 1 mL NaOH
- Observe and record the change in **color** of solution during titration



Clean-Up and Check-Out

- Place the pH electrode in **electrode-cap** filled with **3 M KCl**
- Switch off the pH meter (keep the power cord plugged in)
- Salt solutions resulted from **acid-base neutralization** can be **disposed into the sink**
- Return the **magnetic stir bar** to TA
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- Groups on duty shall stay and help clean up the lab
- This is a **Full Report** experiment:
 - **Have the lab notes and results checked by the TA, and hand in the report next week.**





Data Analysis of Report

1. Average standardized concentration of NaOH

2. Experimental data

3. 3 plots

4. 3 equivalence points

} In Excel

$$C_{\text{NaOH}} \cdot V_{\text{NaOH}} = \frac{\text{Mass}_{\text{KHP}}}{204.22}$$

$$C_{\text{NaOH}} V_{\text{NaOH}} = C_{\text{HOAc}} V_2$$

5. The molar concentration of acetic acid in vinegar (C_{HOAc})

6. Mass percent concentration

- Compare with labels (assume density is same as water)

For example:

$$C_M = 0.737 \text{ M}$$



$$\frac{0.737 \text{ mol/L} \times 60.0 \text{ g/mol}}{1000 \text{ mL/L} \times 1.0 \text{ g/mL}} \times 100\% = 4.4\%$$

7. K_a of acetic acid

- Compare with theoretical value



Example of Data Sheet

Set data to two decimal places

$$V_{\text{titrant}} = V_{\text{read}} - V_i$$

$$0.35 - 0.35 = 0.00$$

$$1.35 - 0.35 = 1.00$$

Mark half- equivalence point

V_{ead} (mL)	V_{titrant} (mL)	pH
0.35	0.00	2.94
1.35	1.00	3.34
2.35	2.00	3.60
3.35	3.00	3.78
4.35	4.00	3.92
5.35	5.00	4.03
6.35	6.00	4.13
7.35	7.00	4.22
8.35	8.00	4.31
9.35	9.00	4.39
10.35	10.00	4.47
11.35	11.00	4.54
12.35	12.00	4.62
13.35	13.00	4.70
14.35	14.00	4.78
15.35	15.00	4.86

Mark equivalence point

V_{ead} (mL)	V_{titrant} (mL)	pH
16.35	16.00	4.95
17.35	17.00	5.05
18.35	18.00	5.17
19.35	19.00	5.30
20.35	20.00	5.47
20.55	20.20	5.51
20.75	20.40	5.56
20.95	20.60	5.60
21.15	20.80	5.66
21.35	21.00	5.72
21.55	21.20	5.79
21.75	21.40	5.85
21.95	21.60	5.93
22.15	21.80	6.04
22.35	22.00	6.17
22.55	22.20	6.33
22.75	22.40	6.57
22.95	22.60	6.86
23.15	22.80	8.35
23.35	23.00	9.20
23.55	23.20	9.80
23.75	23.40	10.21
24.75	24.40	11.09



Example of Data Analysis

$$V_1 = \frac{0.00 + 1.00}{2} = 0.50$$

$$V_2 = \frac{0.50 + 1.50}{2} = 1.00$$

$$\Delta\text{pH}/\Delta V = \frac{3.34 - 2.94}{1.00 - 0.00} = 0.40$$

$$\Delta(\Delta\text{pH}/\Delta V)/\Delta V_1 = \frac{0.26 - 0.40}{1.50 - 0.50} = -0.14$$

Titration curve		First derivative		Second derivative	
V_{NaOH}	pH	V_1	$\Delta\text{pH}/\Delta V$	V_2	$\Delta(\Delta\text{pH}/\Delta V)/\Delta V_1$
0.00	2.94	0.50	0.40	1.00	-0.14
1.00	3.34	1.50	0.26	2.00	-0.08
2.00	3.60	2.50	0.18	3.00	-0.04
.....
22.40	6.57	22.50	1.45	22.60	30.00
22.60	6.86	22.70	7.45	22.80	-16.00
22.80	8.35	22.90	4.25	23.00	-6.25
23.00	9.20	23.10	3.00	23.20	-4.75

Mark
equivalence
point




Example of Titration Curve

Mark

1. equivalence point
2. half- equivalence point


$$V_{\text{eq}} = \frac{19.60 \text{ mL} + 19.40 \text{ mL}}{2} = 19.50 \text{ mL}$$





Example of Figures


Figure 1



Mark half- equivalence point in fig. 1


Mark equivalence point in 3 figures

Figure 2



Mark equivalence point in 3 figures

Figure 3



Mark equivalence point in 3 figures