



# General Chemistry Laboratory

---

## Solubility Product Constant of Silver Acetate



# Preparation


## Collect the following items

- One magnetic stir bar (TA will distribute)
- One 25 mL buret, four 125 mL Erlenmeyer flasks
- One 10 mL graduated pipet and one pipet filler
- Three** glass funnels and one 100 mL beaker
- Ring clamp, Styrofoam cup, NBR gloves



## From your personal equipment

- Four** 100 mL beakers



## Preparation

- Take two clean and dry 100 mL beakers to prepare saturated AgOAc solution
- Wash clean and oven dry two glass funnels ( $\phi$  70) and three 100 mL beakers (label group no. on beaker)



# Notes

**The concentrations of chemicals are changed, which affects the calculation of  $K_{sp}$**

Chemicals (p.134)

- **0.020 M** Potassium thiocyanate, KSCN
- **0.10 M** Silver acetate, AgNO<sub>3</sub>

I. Experimental data (p.137)

- **0.020 M** KSCN



# Objective and Principles


- **Objective:** Determine the solubility product constant ( $K_{sp}$ ) of silver acetate ( $\text{AgOAc}$ ) at specific temperatures
- **Lab techniques:**
  - Gravity filtration
  - Using a graduated pipet
  - Performing titration using a buret
  - Hot plate/magnetic stirrer
  - Lab dispenser






# Solubility Product Constant

- $K_{sp}$  is a specific type of equilibrium constant that describes the dissolution process of salts in water:



$$K_{sp} = [A^{n+}]^m [B^{m-}]^n$$

- When the ion product ( $Q$ ) exceeds  $K_{sp}$ , precipitates would form until the  $Q$  value in the solution equals to  $K_{sp}$






# Solubility Product Constant

- For any specific salt, its  $K_{sp}$  changes with temperature but stays independent of pH and the concentrations of ions in the solution
- Nevertheless, the *solubility* of salts may be affected by the presence of other ions (common ion effect)

## Solubility of AgCl in pure water




Initial		0	0
Change	- x	+ x	+ x
equi.		x	x

$$K_{sp} = [\text{Ag}^+][\text{Cl}^-] = 1.6 \times 10^{-10}$$

$$x^2 = 1.6 \times 10^{-10}$$
$$\rightarrow x = 1.3 \times 10^{-5} \text{ (M)}$$

## Solubility of AgCl in 0.01 M NaCl



Initial		0	<b>0.01</b>
Change	- x	+ x	+ x
equi.		x	<b>0.01 + x</b>




$$K_{sp} = [\text{Ag}^+][\text{Cl}^-] = 1.6 \times 10^{-10}$$

$$x(0.01 + x)$$
$$\simeq (0.01)x = 1.6 \times 10^{-10}$$
$$\rightarrow x = 1.6 \times 10^{-8} \text{ (M)}$$



# Experiment Tasks


1. Prepare saturated silver acetate solution
2. Remove excess AgOAc precipitate
3. Use a standard  $\text{SCN}^-$  solution to titrate and determine  $[\text{Ag}^+]$
4. Based on the known  $[\text{OAc}^-]$ , calculate  $K_{\text{sp}}(\text{AgOAc}) = [\text{Ag}^+][\text{OAc}^-]$
5. Perform the same procedures in both room temperature and a lower temperature





# Step 1: Prepare Saturated AgOAc Solution

- Use a clean and dry 100 mL beaker to take 10.0 mL **0.10 M** of **AgNO<sub>3</sub>(aq)** and 15.0 mL of 0.30 M NaOAc(aq)
- Add the magnetic stirrer into the beaker, stir for **15** minutes
- Record the solution temperature








## Step 2: Filtration to Remove AgOAc Ppt

- Take a 110-mm diameter filter paper, fold it twice and tear off a small piece at one of the two thin outside corners
- Use a ring clamp to set up a glass funnel; expand the filter paper and install it in the funnel
- Perform gravity filtration to remove AgOAc precipitate and collect the filtrate in a 100 mL beaker
- **Rinse** a graduate pipet with few AgOAc(aq) filtrate then transfer assigned portion into 125 mL Erlenmeyer flasks
  - **5.00 mL** (coarse titration)
  - **10.00 mL** (fine titration)
- Add **1 mL of Fe<sup>3+</sup> indicator**







# Step 3: Titration of $[Ag^+]$



- Use a clean and dry 100 mL beaker to take **40.0 mL 0.020 M KSCN(aq)**
- Rinse the buret twice with KSCN solution, then fill the buret
- Titrate the AgOAc solution with **0.020 M KSCN** until the solution appears a light orange color (read and record  $V_i$  and  $V_f$  to 0.01 mL)





Over  
titration

Correct  
end point




## Step 4: $K_{sp}$ at Low Temperature

- Use another clean and dry 100 mL beaker and take 10.0 mL **0.10 M** of  $\text{AgNO}_3(\text{aq})$  and 15.0 mL of 0.30 M  $\text{NaOAc}(\text{aq})$
- Stir and mix the solution in an ice-water bath for **15** minutes, record the temperature of solution
- Repeat Steps 2 & 3 to determine  $[\text{Ag}^+]$  in saturated  $\text{AgOAc}(\text{aq})$





# Calculate $[OAc^-]$



Before mixing     $10.0 \text{ mL} \times 0.10 \text{ M} = 1.0 \text{ mmol}$      $15.0 \text{ mL} \times 0.30 \text{ M} = 4.5 \text{ mmol (excess)}$



After mixing     $1.0 - 1.0 = 0 \text{ mmol}$      $4.5 - 1.0 = 3.5 \text{ mmol}$      $1.0 \text{ mmol}$   
Ppt. occurs

- Because the  $K_{sp}$  of AgOAc is quite small, the precipitation reaction is assumed to be completed with NaOAc being the excess reagent
- Total volume after mixing = 25.0 mL

$$\text{Excess reagent: } [OAc^-] = \frac{3.5 \text{ mmol}}{25.0 \text{ mL}} = 0.14 \text{ M}$$



# Calculate $[Ag^+]$ (Volhard Method)



- $SCN^-$  reacts with  $Ag^+$  first
- At titration end point, # moles  $Ag^+$  = # moles  $SCN^-$
- $[Ag^+] \times V_1 = [SCN^-] \times (V_f - V_i)$



# Calculate $K_{sp}$ of AgOAc



Concentration before  
Equilibrium

0 M

0.14 M

Concentration after  
**Dissolution** equilibrium

**x** M

(0.14 + **x**) M

$$K_{sp} = [\text{Ag}^+][\text{OAc}^-] = \mathbf{x (0.14 + x)}$$



# Additional Notes


- Use only clean and dry beakers for mixing  $\text{AgNO}_3$  and  $\text{NaOAc}$
- Rinse the graduated pipet with sat.  $\text{AgOAc}(\text{aq})$ , and rinse the buret with  $\text{KSCN}(\text{aq})$  to ensure that  $[\text{Ag}^+]$  can be accurately determined
- Take only **40 mL 0.020 M KSCN(aq)** by lab dispenser
- Add **1 mL  $\text{Fe}^{3+}$**  indicator for each titration (not 1 drop)
- Temperature of  $\text{AgOAc}(\text{aq})$  should be recorded right before the gravity filtration step
- Avoid contacting  $\text{AgNO}_3(\text{aq})$  or  $\text{AgOAc}(\text{aq})$  with skin (wear NBR gloves at all times)
- **Be careful in not losing the magnetic stir bar into the sink or waste bin**



# Clean-Up and Check-Out

## Brief Report

- Ag-containing solution and precipitates should be disposed into the designated heavy metal waste bin
- The cleaned buret and transfer pipet should be clamped upside-down on the buret clamp
- Wash Erlenmeyer flasks thoroughly with brush and detergent
- Return the magnetic stir bar to TA
- Clean up the lab bench and check personal equipment inventory (have an associate TA sign the check list)
- This is a **Brief Report** experiment:
- Groups on duty shall stay and help clean up the lab







# Example of Experimental Data

Test item		Room temp. (°C) (ex. 24.0 °C)		Ice bath (°C) (ex. 3.5 °C)	
		Coarse titration	Fine titration	Coarse titration	Fine titration
Volume of saturated AgOAc (mL)		5.00	10.00	5.00	10.00
0.020 M KSCN	$V_i$ (mL)				
	$V_f$ (mL)				
Titration volume	$\Delta V$ (mL)				




# Calculate $[Ag^+]$ , $[OAc^-]$ , and $K_{sp}$

Temp. ( $^{\circ}C$ )	(1) 24.0	(2) 3.5
$[Ag^+]$ (M)	$[Ag^+] \times 10.00 = 0.020 \times 11.11$ $[Ag^+] = 0.023$	
$[OAc^-]$ (M)	$0.14 + 0.023 = 0.16$	
$K_{sp}$	$3.7 \times 10^{-3}$	

\* Compare the  $K_{sp}$  of silver acetate at different temperature



# T2 – Stirrer/Hot Plate



- Connect the stirrer/hot plate to a grounded 110 V power outlet (replace damaged power cord and plug immediately)
- Keep power cord away from the ceramic heating top
- Clean the heating top with non-corrosive detergent after use or when liquid spills
- NEVER heat a large amount of volatile and flammable liquid (e.g. ether, acetone) directly on the hot plate
- If the stirring bar jumps erratically, turn the stirring function off and adjust the vessel position, then restart the stirring
- Do not remove the stirring bar from solution with hand – instead use a Teflon-coated magnetic rod (“fishing pole”)



# T6 – Gravity Filtration




Figure T6-1 Folding a filter cone

- Fold a round filter paper in half for two times. Tear off a small piece at one of the two thin outside corners
- Expand the filter paper (from the intact fold) into a cone shape. Fit the filter paper into a funnel (the edge of filter paper should not exceed the top of funnel)
- Use a ring clamp to support the glass funnel. The tip of the funnel should touch the sidewall of the receiving vessel
- Pour the liquid into the paper cone (not on the glass funnel). Use a glass rod to decant the liquid
- Fill the paper cone no more than 2/3 full
- After filtration, use a tweezer to separate the filter paper from funnel (don't use hand)





Figure T6-2 Setup of gravity filtration



# T12.2 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 1

- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid to the 0.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Press valve **E** to drain the liquid to the 5.00 mL marking
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use






# T12.3 – Measuring (Graduated) Pipet

Deliver 5.00 mL solution – Method 2




- Clean a 10 mL pipet and rinse it twice with small amount of the liquid to be transferred
- Press valve **A** of the pipet filler and simultaneously squeeze the bulb to expel air from it, then insert the top of pipet gently into the pipet filler
- Bring the pipet tip below the liquid surface, press valve **S** to draw liquid until it rises above the 5.00 mL marking
- Remove the pipet filler and quickly use an index finger to close the top of pipet. Use the finger to adjust the liquid level to the 5.00 mL marking
- Wipe off any excess liquid near the pipet tip
- Use the other hand to hold the new container. Maintain the pipet in a vertical position and let its tip touch the inner wall of the container. Release the index finger so that liquid is transferred
- Do not force out any liquid remaining at the tip
- Wash the pipet thoroughly after use





# T14 – Titration

- Clean the buret with DI water, then rinse twice with ~5 mL of titrant (use a funnel to add titrant)
- Open the stopcock to repel the air at the buret tip
- Adjust the height of buret so that its tip is lower than the lid of receiving flask
- Read and record the initial volume ( $V_i$ ) on the buret to 0.01 mL
- With the stopcock on the right side, use your left hand to control the stopcock while the right hand swirls the receiving flask in a circular motion
- At the titration end point, read and record the final volume ( $V_f$ ) on the buret to 0.01 mL
- After the experiment, wash the buret and let it dry upside-down on the buret clamp





# Lab Dispenser

- 1) Check the pre-set volume on the dispenser. Do not change the setting unless instructed to do so
- 2) Place the receiving flask under the tip of dispenser
- 3) To remove the air bubbles in the dispenser, lightly pull the piston pump up and down for several times
- 4) Gently pull the piston pump up until it reaches the end of travel range, then slowly push the piston down to obtain the solution

