



General Chemistry Laboratory

Recrystallization & Melting Point Determination



Preparation

Collect the following items

- ☐ Melt-Temp (analog) or MP-2D (digital)
- ☐ Nine capillary tubes
- ☐ Hollow glass tube
- ☐ Büchner funnel
- ☐ Filter paper
- ☐ Your benzoic acid and acetanilide products

From your personal equipment

- ☐ Capillary tube stand
- ☐ Hot plate
- ☐ Two 50 mL Erlenmeyer flasks
- ☐ Suction flask
- ☐ Rubber adapter cone
- ☐ Water aspirator



Capillary tube stand



Objective and Principles

- **Objective:** Use recrystallization to purify compound and determine the melting point (mp)
- **Lab techniques**
 - Recrystallization
 - Suction filtration
 - Packing capillary tubes
 - Measuring melting points of organic compounds
- **Tasks**
 1. Recrystallization of the crude benzoic acid
 2. Measuring the mp of reference compounds
 3. Measuring the mp of crude and purified compounds



Principle

■ Melting point

- The temperature at which substance changes from solid to liquid at atmospheric pressure
- When a pure crystalline substance melts, the melting point range should not exceed 1°C . An impure substance shows a larger melting point range than a pure substance
- The melting point can be used to
 - ✓ determine the identity of a known compound
 - ✓ define an unknown compound at later times
 - ✓ determine the purity of a substance



Principle

■ Recrystallization

- Dissolve the sample in a suitable hot solvent. After the solution is cooled, a supersaturated solution will be produced, and the solute will crystallize out
- During recrystallization, the sample adheres to the lattice in an orderly manner, without any impurities included that increases the purity
- Choose suitable solvent
 - has high solubility for the sample at high temperature, but low solubility at low temperature
 - has high solubility for the impurity; hence, it will not be crystallized out
 - non-flammable, non-toxic, cheap, and volatile...



Step 1: Packing Capillary Tube



Sealed end down and open end up

- Pack 2 capillary tubes for each sample:
 - (1) benzoic acid (mp 122°C)
 - (2) acetanilide (mp 113°C)
 - (3) benzoic acid/acetanilide mixture (1:1)

- Transfer a dry and finely powdered sample onto a piece of weighing paper
- Insert the open-end of a capillary tube into the pile of the sample
- Knock the closed-end of the capillary tube on the bench top 3-4 times
- The packed sample is *ca.* 2 mm height



Step 2: Determine the Melting Point Range



Heating rate:

Initial: 10~15°C/min

Lower than mp 10~20°C: 2°C/min


Lower than mp 2~4°C: 1°C/min

✓ Line up the samples from the one with lowest mp and end with the highest one

- Set the voltage control knob to zero, then turn on the power
- To pre-measure a rough mp range for three samples, apply a faster heating rate of about 5°C/min
- Cool the Melt-Temp apparatus to lower than the approximate mp about 15°C
- Replace the capillary tubes and slowly heat the samples at a rate of 2°C/min to determine the mp
- Record the mp range of the sample, i.e. the temperature range from starting to melt to completely melted



Step 3: Recrystallization





- Use a 50 mL Erlenmeyer flask to heat and boil some DI water
- Prepare 2 capillary tubes of crude benzoic acid for mp determination
- Transfer the remaining sample to a 50 mL Erlenmeyer flask after weighing
- Use dropper to add some hot DI water to the flask
- Heat and swirl the solution gently to dissolve sample
- Cool it down slowly to room temp., and let it stand to grow crystals

- ✓ Calculate the amount of hot water needed base on solubility
- ✓ As water is the solvent, heat the solution on hot plate directly instead of water bath



Step 4: Collect Product and Determine the Melting Point



- Use an ice water bath to cool the mixture to a lower temperature and produce more crystals
- Use suction filtration to collect crystals
- Wash the crystals with small portion of cold-water
- Suction dry for 10 min.
- Collect the crystals on filter paper and press to dry
- Weigh the dried crystals and determine the yield
- Determine mp of (1) crude benzoic acid, (2) recrystallized benzoic acid, (3) crude acetanilide



Additional Notice

Melting point determination

- The sample inside the capillary tube should be 2 mm in height
- Glass tube for knocking the capillary tube should be kept clean and dry; wash and oven dry after class
- Record the melting point range of sample
- Replace the capillary tubes in the second run measurement

Recrystallization

- Use an Erlenmeyer flask for recrystallization
- Do not use flame or hot plate to directly heat flammable organic solvents. Use water bath instead
- Leave the solution to cool to room temperature, and let it stand to grow crystals slowly that increases the purity



Clean-Up and Check-Out



- Dispose of benzoic acid, acetanilide, and capillary tubes to designated waste bins
- Clean up the lab bench and check personal equipment inventory (have an associate TA signed the check list)
- This is a **Brief Report** experiment:
 - Hand in prelab/lab note/report together to the TA
- Groups on duty shall stay and help clean up the lab